

Technical specifications Eagle 620 - 625

Eagle 620

total height	total length	total width (closed)	total width (opened)	height arm	users weight	total weight
174 cm	111,5 cm	55 cm	103,5 cm	from 58 to 160cm	150 kg	37 kg

Deviations in sizes and weight (1,5 cm and 1,5 kg) are possible

Eagle 625

total height	total length	total width (closed)	total width (opened)	height arm	users weight	total weight
184 cm	129 cm	68 cm	122,5 cm	from 65 to 173cm	175 kg	42 kg

Deviations in sizes and weight (1,5 cm and 1,5 kg) are possible

Your guarantee of quality.

Regarding quality, the **Patientlift** complies with the highest norms, laid down by certified test institutes in most European countries.

N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be



VERMEIREN
Opening up new horizons in mobility

For more information, contact your local dealer



R.E.: N.V. Vermeiren N.V., Vermeirenplein 1/15 - 2920 Kalmthout - Belgium - 08/2010

Subject to technical alterations

EAGLE

620 - 625



PATIENTLIFT

The Vermeiren patientlift was designed to bring a real comfort to the patient during the transfers but also for the people who are taking care of the patient.



Optional:
U-shaped slings
M for model 620
L for model 625

The patientlift Eagle 620 and Eagle 625 are equipped with the Linak HB7X engine feeding from a battery with an autonomy of approximately 45 lifts.



Security and simplicity...



Fast and easy folding system, effortlessly without any tool thanks to the wing screws and its axis with fast disassembling.

Easy transport and storage



Perfect stability is always ensured when you want to adjust the space between the two legs of the base by means of a very simple system of pedals.



The low height of the base allows it to slide under most beds on the market.

Manual emergency lowering model Eagle 620



Manual emergency lowering model Eagle 625



Remote control Linak HB7X (IPX6)



Security button



Electric emergency lowering of the arm




Vermeiren designed this patientlift with specific attention to functionality, safety and reliability. Quality wheels of 100 mm ensure that this patientlift can be used on almost any underground.

SLING with HEADSUPPORT

available in
M - L - XL (solid & net)


	A	B	C	D	E
M	78	40	45	30	39
L	84	45	50	36	41
XL	91	51	56	44	42



SLING with TOILETOPENING

available in
M - L - XL (solid & net)

	A	B	C	D
M	78	25	60	30
L	84	26	70	36
XL	91	27	82	44



STANDARD SLING

available in
M - L - XL (solid & net)

	A	B	C	D
M	78	40	45	30
L	84	45	50	36
XL	91	51	56	44

